

Home Care Series

COMPRESSOR NEBULIZER HC - M3

- ▶ Short inhalation time, fast relief – with extensive accessories
- ▶ For treating the upper and lower airways – with compressed-air technology
- ▶ High proportion of respirable particles–helps with colds, asthma and respiratory diseases

Medical
nebulizer



Low
noise



Small atomized
particles



Safe and
Durable



COMPRESSOR NEBULIZER HC - M2

PRODUCT Features

- Lightweight and portable, atomize anytime, anywhere.
 - 60dB soft sound atomization, multiple noise reduction, soft sound does not disturb.
 - Adjustable atomization cup, the amount of fog can be adjusted as you like.
 - The atomized particles are fine and ultra-fine atomized, reaching deep into the lungs.
 - Mask/mouthpiece, two modes selectable.
 - Medical grade pp material, can directly contact with the baby's skin.
 - 220V low voltage DC, safety protection.
- The atomization rate and spray volume are constant, gentle and not throat-choking.

PRODUCT SPECIFICATION

Power supply	AC 220V
Atomization rate	≥0.2mL/min
Residual liquid amount	≤1.0 mL
Median particle size	2.5um±25%
Operating noise level	<45dB
Equivalent volume size distribution	
of fog particles	≥70%
Service life	3 years
Working power	<100W
Power supply	2AA batteries or DC5V
Net weight	About 1250g (without battery)
Dimension	150*150*90mm



Home Care Series



HANDHELD NEBULIZER HC-N2

PRODUCT INFORMATION

Atomizing rate	0.35ml/min
Median particle size	3um±25%
Medicine cup capacity	8ml
Operating noise level	<50dB(A)
Nominal frequency	108kHz
Residual liquid amount	≤0.5ml
Power supply	2"AA batteries or DC5V (optional lithium battery)
Net weight	about 109.5g (no battery)
Dimension	130*48*50mm

PRODUCT Features

- Compact and portable, low noise design.
- Atomized particles can reach deep into the lungs.
- High sensitivity water detection, small liquid residue is not wasted.
- Intelligent low power mode with long standby time.
- Self-contained cleaning mode.



ARM TYPE BLOOD PRESSURE MONITOR HC-BPA2

PRODUCT Features

Blood pressure measure unit:

mmHg/kPa. According to habit to switch.

Hierarchical display blood pressure measurement result.

Low battery warning.

Overpressure protection, automatic fast exhaust when air pressure over 295mmHg(20ms).

Automatic power off when there is no operation for 2 minutes.

Streamlined appearance design.

Accommodates two users, storing up to 99 readings per user on one device.

Optional function: Voice, backlight, Bluetooth, lithium battery.

PRODUCT SPECIFICATION

NIBP	Static pressure measurement range: 0-280 mmHg
PR	Accuracy:±3mm Hg(±0.4kPa) Measurement range:40BPM~199BPM
Measurement	upper arm
Display	segment LCD display
Power	4 AA alkaline batteries, DC6.0V
Product dimension	130.9*96.1*63.1mm
Product weight	237g



Home Care Series

DAILY OXYGEN PROTECTION AND HEALTH CARE

Model:HC-OXY502N



Imported molecular sieve gives you a full sense of
security with natural enjoyment!

homecare —Homecare Oxygen Concentrator, protect the family's oxygen health.